

Product Datasheet

Rat LTF ELISA Kit (orb440638)

Description	Rat Lactoferrin (LTF) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D3ZAB1
Sensitivity	12.7pg/mL
Entrez	301034
Expiration Date	Please enquire.