

## Product Datasheet

### Human MAdCAM1 ELISA Kit (orb440648)

<b>Description</b>	Human Mucosal Addressin Cell Adhesion Molecule 1 (MAdCAM1) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, urine and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, urine and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q13477</a>
<b>Sensitivity</b>	0.125ng/mL
<b>Entrez</b>	<a href="#">8174</a>
<b>Expiration Date</b>	Please enquire.