

Product Datasheet

Human MAP1LC3a ELISA Kit (orb440659)

Description	Human Microtubule Associated Protein 1 Light Chain 3 Alpha (MAP1LC3a) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9H492
Sensitivity	0.063ng/mL
Entrez	84557
Expiration Date	Please enquire.