

Product Datasheet

Rat MAPt ELISA Kit (orb440673)

Description	Rat Microtubule Associated Protein Tau (MAPt) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P19332
Sensitivity	25pg/mL
Entrez	29477
Expiration Date	Please enquire.