

## Product Datasheet

### Human MASP1 ELISA Kit (orb440676)

<b>Description</b>	Human Mannose Associated Serine Protease 1 (MASP1) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	93.75-6000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma and other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<b>P48740</b>
<b>Sensitivity</b>	32.12pg/mL
<b>Entrez</b>	<b>5648</b>
<b>Expiration Date</b>	Please enquire.