

Product Datasheet

Mouse MASP1 ELISA Kit (orb440677)

Description	Mouse Mannose Associated Serine Protease 1 (MASP1) Sandwich ELISA Kit
Reactivity	Mouse
Range	2.343-150ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P98064
Sensitivity	0.8ng/mL
Entrez	17174
Expiration Date	Please enquire.