

Product Datasheet

Monkey MCP1 ELISA Kit (orb440710)

Description	Monkey Monocyte Chemotactic Protein 1 (MCP1) Sandwich ELISA Kit
Reactivity	Monkey
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P61274
Sensitivity	0.054ng/mL
Entrez	102135739
Expiration Date	Please enquire.