


Product Datasheet

CD25 antibody (Biotin) (orb44072)

Description	Mouse monoclonal antibody against CD25 conjugated to Biotin
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, IP
Immunogen	PHA-activated peripheral blood leucocytes
Target	CD25
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P01589
Entrez	3559
Dilution Range	Flow cytometry: Recommended dilution: 1 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (surface stainin...