

Product Datasheet

Mouse MDC ELISA Kit (orb440726)

Description	Mouse Macrophage Derived Chemokine (MDC) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9R1V4
Sensitivity	5.8pg/mL
Entrez	11488
Expiration Date	Please enquire.