

## Product Datasheet

### Chicken MIP1b ELISA Kit (orb440784)

<b>Description</b>	Chicken Macrophage Inflammatory Protein 1 Beta (MIP1b) Sandwich ELISA Kit
<b>Reactivity</b>	Gallus
<b>Range</b>	62.5-4000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>Q90826</b>
<b>Sensitivity</b>	25.8pg/mL
<b>Entrez</b>	<b>395551</b>
<b>Expiration Date</b>	Please enquire.