

Product Datasheet

Monkey MIP1b ELISA Kit (orb440790)

Description	Monkey Macrophage Inflammatory Protein 1 Beta (MIP1b) Sandwich ELISA Kit
Reactivity	Monkey
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8HYQ2
Sensitivity	9.5pg/mL
Entrez	574217
Expiration Date	Please enquire.