

Product Datasheet

Mouse MIP3a ELISA Kit (orb440794)

Description	Mouse Macrophage Inflammatory Protein 3 Alpha (MIP3a) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O89093
Sensitivity	27pg/mL
Entrez	20297
Expiration Date	Please enquire.