

Product Datasheet

Canine MIP3b ELISA Kit (orb440796)

| | |
|--------------------------|--|
| Description | Canine Macrophage Inflammatory Protein 3 Beta (MIP3b) Sandwich ELISA Kit |
| Reactivity | Canine |
| Range | 15.625-1000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma and other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | Q68AY7 |
| Sensitivity | 6.3pg/mL |
| Entrez | 448793 |
| Expiration Date | Please enquire. |