

Product Datasheet

Mouse MMP15 ELISA Kit (orb440832)

Description	Mouse Matrix Metalloproteinase 15 (MMP15) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O54732
Sensitivity	28.4pg/mL
Entrez	17388
Expiration Date	Please enquire.