

Product Datasheet

Canine MMP3 ELISA Kit (orb440859)

Description	Canine Matrix Metalloproteinase 3 (MMP3) Sandwich ELISA Kit
Reactivity	Canine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q6Y4Q5
Sensitivity	0.063ng/mL
Entrez	403445
Expiration Date	Please enquire.