

Product Datasheet

Rat MMRN1 ELISA Kit (orb440885)

Description	Rat Multimerin 1 (MMRN1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D4A3E0
Sensitivity	0.113ng/mL
Entrez	500152
Expiration Date	Please enquire.