

Product Datasheet

Human MRP1 ELISA Kit (orb440915)

Description	Human Motility Related Protein (MRP1) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P33527
Sensitivity	12.4pg/mL
Entrez	4363
Expiration Date	Please enquire.