

Product Datasheet

Human MX1 ELISA Kit (orb440984)

Description	Human Myxovirus Resistance 1 (MX1) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P20591
Sensitivity	0.27ng/mL
Entrez	4599
Expiration Date	Please enquire.