

Product Datasheet

Human MYH3 ELISA Kit (orb441002)

Description	Human Myosin Heavy Chain 3, Skeletal Muscle, Embryonic (MYH3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P11055
Sensitivity	0.126ng/mL
Entrez	4621
Expiration Date	Please enquire.