

## Product Datasheet

### Human NEF3 ELISA Kit (orb441072)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Neurofilament 3 (NEF3) Sandwich ELISA Kit                                          |
| <b>Reactivity</b>        | Human                                                                                    |
| <b>Range</b>             | 15.625-1000pg/mL                                                                         |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | serum, plasma, tissue homogenates, cerebrospinal fluid and other biological fluids.      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cerebrospinal fluid and other biological fluids.      |
| <b>Uniprot ID</b>        | <b>P07197</b>                                                                            |
| <b>Sensitivity</b>       | 6.2pg/mL                                                                                 |
| <b>Entrez</b>            | <b>4741</b>                                                                              |
| <b>Expiration Date</b>   | Please enquire.                                                                          |