

## Product Datasheet

### Human NEFL ELISA Kit (orb441075)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Human Neurofilament, Light Polypeptide (NEFL) Sandwich ELISA Kit                         |
| <b>Reactivity</b>        | Human  |
| <b>Range</b>             | 15.625-1000pg/mL   |
| <b>Form/Appearance</b>   | Traditional  |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | serum, plasma, cerebrospinal fluid and other biological fluids.                          |
| <b>Sample Types</b>      | serum, plasma, cerebrospinal fluid and other biological fluids.                          |
| <b>Uniprot ID</b>        | <a href="#">P07196</a>   |
| <b>Sensitivity</b>       | 5.9pg/mL   |
| <b>Entrez</b>            | <a href="#">4747</a>   |
| <b>Expiration Date</b>   | Please enquire.  |