

Product Datasheet

General Noradrenaline ELISA Kit (orb441077)

Description General Noradrenaline (NE) competitive inhibition ELISA Kit

Reactivity All

Range 78.125-5000pg/mL

Form/Appearance Traditional

Storage All reagents should be kept according to the instructions on the packaging upon receipt.

Note For research use only

Application notes serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.

Sample Types serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.

Sensitivity 27.2pg/mL

Expiration Date Please enquire.