

Product Datasheet

Mouse NMB ELISA Kit (orb441140)

Description	Mouse Neuromedin B (NMB) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D3Z621
Sensitivity	4.75pg/mL
Entrez	68039
Expiration Date	Please enquire.