

## Product Datasheet

### Canine NOS2 ELISA Kit (orb441156)

<b>Description</b>	Canine Nitric Oxide Synthase 2, Inducible (NOS2) Sandwich ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	0.312-20ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">O62699</a>
<b>Sensitivity</b>	0.126ng/mL
<b>Entrez</b>	<a href="#">403822</a>
<b>Expiration Date</b>	Please enquire.