

Product Datasheet

Human NOTCH2NL ELISA Kit (orb441169)

Description	Human Notch2 N-Terminal Like Protein (NOTCH2NL) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q7Z3S9
Sensitivity	0.117ng/mL
Entrez	388677
Expiration Date	Please enquire.