

Product Datasheet

Rat NPS ELISA Kit (orb441192)

Description	Rat Neuropeptide S (NPS) competitive inhibition ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P0C0P7
Sensitivity	4.75pg/mL
Entrez	100360071
Expiration Date	Please enquire.