

Product Datasheet

Rat NPY1R ELISA Kit (orb441203)

| | |
|--------------------------|--|
| Description | Rat Neuropeptide Y Receptor Y1 (NPY1R) Sandwich ELISA Kit |
| Reactivity | Rat |
| Range | 0.312-20ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | tissue homogenates and other biological fluids. |
| Sample Types | tissue homogenates and other biological fluids. |
| Uniprot ID | P21555 |
| Sensitivity | 0.128ng/mL |
| Entrez | 29358 |
| Expiration Date | Please enquire. |