

Product Datasheet

Rat NPY2R ELISA Kit (orb441206)

Description	Rat Neuropeptide Y Receptor Y2 (NPY2R) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9ERCO
Sensitivity	0.054ng/mL
Entrez	66024
Expiration Date	Please enquire.