

Product Datasheet

Mouse NPY5R ELISA Kit (orb441208)

Description	Mouse Neuropeptide Y Receptor Y5 (NPY5R) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	070342
Sensitivity	0.057ng/mL
Entrez	18168
Expiration Date	Please enquire.