

## Product Datasheet

### Mouse NT-ProANP ELISA Kit (orb441261)

<b>Description</b>	Mouse N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) Sandwich ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	78.125-5000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Sensitivity</b>	29pg/mL
<b>Expiration Date</b>	Please enquire.