

Product Datasheet

Porcine OPN ELISA Kit (orb441335)

Description	Porcine Osteopontin (OPN) Sandwich ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P14287
Sensitivity	29.8pg/mL
Entrez	397087
Expiration Date	Please enquire.