

Product Datasheet

Human P2RY12 ELISA Kit (orb441370)

Description	Human Purinergic Receptor P2Y, G Protein Coupled 12 (P2RY12) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9H244
Sensitivity	0.063ng/mL
Entrez	64805
Expiration Date	Please enquire.