

Product Datasheet

Bovine PGA ELISA Kit (orb441523)

Description	Bovine Pepsinogen A (PGA) Sandwich ELISA Kit
Reactivity	Bovine
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P00792
Sensitivity	0.58ng/mL
Entrez	414350
Expiration Date	Please enquire.