

Product Datasheet

Rat PGA ELISA Kit (orb441526)

Description	Rat Pepsinogen A (PGA) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9JJX2
Sensitivity	15.5pg/mL
Entrez	60372
Expiration Date	Please enquire.