

Product Datasheet

Mouse PGR ELISA Kit (orb441544)

Description	Mouse Progesterone Receptor (PGR) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	A6H6A5
Sensitivity	0.128ng/mL
Entrez	18667
Expiration Date	Please enquire.