

Product Datasheet

Human PIGR ELISA Kit (orb441559)

Description	Human Polymeric Immunoglobulin Receptor (PIGR) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine, saliva and other biological fluids.
Sample Types	serum, plasma, urine, saliva and other biological fluids.
Uniprot ID	P01833
Sensitivity	1.25ng/mL
Entrez	5284
Expiration Date	Please enquire.