

Product Datasheet

Human PIK3C2a ELISA Kit (orb441573)

Description	Human Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide (PIK3C2a) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200U/L
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	000443
Sensitivity	1.6U/L
Entrez	5286
Expiration Date	Please enquire.