

Product Datasheet

Human PIK3Cb ELISA Kit (orb441575)

Description	Human Phosphoinositide-3-Kinase Catalytic Beta Polypeptide (PIK3Cb) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P42338
Sensitivity	0.059ng/mL
Entrez	5291
Expiration Date	Please enquire.