

Product Datasheet

Human PKM2 ELISA Kit (orb441606)

Description	Human Pyruvate Kinase, Muscle (PKM2) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P14618
Sensitivity	29pg/mL
Entrez	5315
Expiration Date	Please enquire.