

Product Datasheet

Rat PLA2G2A ELISA Kit (orb441624)

Description	Rat Phospholipase A2, Group IIA (PLA2G2A) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P14423
Sensitivity	0.116ng/mL
Entrez	29692
Expiration Date	Please enquire.