

Product Datasheet

Mouse PMP22 ELISA Kit (orb441669)

Description	Mouse Peripheral Myelin Protein 22 (PMP22) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P16646
Sensitivity	0.058ng/mL
Entrez	18858
Expiration Date	Please enquire.