

## Product Datasheet

### Human POTE<sup>G</sup> ELISA Kit (orb441690)

<b>Description</b>	Human POTE Ankyrin Domain Family, Member G (POTE <sup>G</sup> ) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Sample Types</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q6S5H5</a>
<b>Sensitivity</b>	0.117ng/mL
<b>Entrez</b>	<a href="#">404785</a>
<b>Expiration Date</b>	Please enquire.