

Product Datasheet

Rat PPAR γ C1a ELISA Kit (orb441700)

Description	Rat Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha (PPAR γ C1a) Sandwich ELISA K...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9QYK2
Sensitivity	0.061ng/mL
Entrez	83516
Expiration Date	Please enquire.