

Product Datasheet

Mouse Pancreatic Polypeptide ELISA Kit (orb441718)

Description	Mouse Pancreatic Polypeptide (PP) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P10601
Sensitivity	8.55pg/mL
Entrez	19064
Expiration Date	Please enquire.