

Product Datasheet

Human PRDX4 ELISA Kit (orb441731)

Description	Human Peroxiredoxin 4 (PRDX4) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q13162
Sensitivity	12.3pg/mL
Entrez	10549
Expiration Date	Please enquire.