

Product Datasheet

Mouse PROZ ELISA Kit (orb441771)

Description	Mouse Protein Z (PROZ) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q9CQW3
Sensitivity	6.5pg/mL
Entrez	66901
Expiration Date	Please enquire.