

Product Datasheet

Rat PTGES ELISA Kit (orb441810)

Description	Rat Prostaglandin E Synthase, Microsomal (PTGES) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9JHF3
Sensitivity	0.068ng/mL
Entrez	59103
Expiration Date	Please enquire.