

Product Datasheet

Rat PTGS2 ELISA Kit (orb441817)

Description	Rat Prostaglandin Endoperoxide Synthase 2 (PTGS2) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35355
Sensitivity	25.7pg/mL
Entrez	29527
Expiration Date	Please enquire.