

Product Datasheet

Human PTH ELISA Kit (orb441820)

Description	Human Parathyroid Hormone (PTH) competitive inhibition ELISA Kit
Reactivity	Human
Range	1.56-100pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01270
Sensitivity	0.49pg/mL
Entrez	5741
Expiration Date	Please enquire.