

Product Datasheet

Mouse PTMS ELISA Kit (orb441832)

Description	Mouse Parathymosin (PTMS) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	46.875-3000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9D0J8
Sensitivity	14.2ng/mL
Entrez	69202
Expiration Date	Please enquire.